


NEW ACE ANTERIOR / DIENTES ANTERIORES NEW ACE

UPPER / SUPERIOR

		A	B	C
▲				
T1		42.6	7.6	10.4
▲				
T2		41.9	7.6	9.9
▲				
T3		45.1	8.2	11.4
▲				
T4		45.0	8.1	10.7
▲				
T5		46.5	8.4	11.3
▲				
T6		48.2	8.6	11.6
▲				
TL4		44.7	8.5	10.9
▲				
TL5		46.9	8.8	11.1
▲				
TL6		49.4	9.1	11.7
▲				
TL7		51.7	9.5	12.3
■				
S2		43.0	7.7	9.5
■				
S3		43.9	8.1	9.8

■				
S4		45.1	8.0	10.6
■				
S5		47.5	8.5	10.5
■				
S6		47.7	8.8	11.4
■				
S7		49.9	9.3	12.6
■				
S8		52.3	9.8	12.9
■				
SS2		41.9	7.8	9.1
■				
SS3		44.3	8.0	9.0
●				
O2		41.4	7.3	9.3
●				
O3		44.1	8.0	10.0
●				
O4		46.6	8.4	10.0
●				
O5		48.8	8.7	10.4



A B C